



EXTRACTION PROTOCOL


© 2002 L. Saphonia Gee

1. Cleanse *Apply pre-cleanser solution for one minute to soften the sebum (oil) and to decrease the incidence of blockage within the follicle openings.*

2. Steam *Apply a warm compress to the targeted area or shower for one minute to open the pores.*

3. Extraction Pads *Wrap fingers with pads*

4. Treatment *First stretch the skin apart and then mound the skin in the targeted area. Gently manipulate with alternating side to side and top to bottom movements extracting the oil.*



5. Solution *Apply the post-extraction solution to close the pores and reduce the incidence of bacteria from entering the post-extraction site.*

6. Post-Extraction Mask *Apply and let remain for 5 to 10 minutes. Mask will re-hydrate to calm the site, recondition the skin to reduce inflammation, and even out skin texture.*

7. Remove the Mask *Apply cold water,*

8. Reapply *both post-extraction solution followed by your daily spf lotion.*

Notes:
